
Fusion mixture, Hi-LR[®]

RM10821

Product Identifier

HS Code : 3822 90 90
Storage : Below 30°C
Shelf life : 4 years

Technical Specification

Appearance : White to off-white hygroscopic crystals or powder or solid
Solubility : 33.3 mg soluble in 1 mL of water

Safety Information

Hazard Pictogram(s) :



Signal Word : Warning
Hazard Statement(s) : H302- H315- H319- H335
Precautionary Statement(s) : P261- P280- P312- P403 + P233- P501
UN No. : Not dangerous goods
Class : -
Packing Group : -
RTECS : Not available
WGK : 1