

RM10823 Span 80®**Product Number Packing**

RM10823	: 2.5L
RM10823	: 500ML

Product Information

Product Code	: RM10823
Product Name	: Span 80®
Synonym	: Sorbitan monooleate
Molecular Formula	: C ₂₄ H ₄₄ O ₆
Molecular Weight	: 428.62
CAS No.	: 1338-43-8
EC No.	: 215-665-4
HS Code	: 3824 90 97
Other Information	: Nonionic surfactant
Shelf Life	: 4 years

Technical Specification

Appearance	: Yellow to brownish yellow viscous liquid
Solubility	: 1 mL miscible in 1 mL of ethanol
Acid value (mg KOH/mg)	: <= 8.00
Saponification value (mg KOH/mg)	: 145 - 160
Hydroxyl value (mg KOH/mg)	: 190 - 210

Risk and Safety Information

WGK	: 1
RTECS	: WG2932400
Flash Point(°F)	: >235.4 °F
Flash Point(°C)	: >113 °C
Storage Temperature(°C)	: Store below 30°C

Transport Information

Marine Pollutant	: No
ADR/RID	: Not Dangerous Goods
IMDG	: Not Dangerous Goods
IATA	: Not Dangerous Goods

Key : ® : Span is a registered trademark of Croda International PLC