

RM10825 Linseed oil, L.R.**Product Number Packing**

RM10825 : 500ML

Product Information

Product Code : RM10825
Product Name : Linseed oil, L.R.
CAS No. : 8001-26-1
EC No. : 232-278-6
HS Code : 1515 19 90
Shelf Life : 4 years

Technical Specification

Appearance : Yellow to brownish-yellow clear liquid (turns dark and gradually thickens on air exposure)
Solubility : 0.1 mL miscible in 5 mL of chloroform
Refractive index (n_{20/D}) : 1.4750 - 1.4850
Density (at 25°C) : 0.926 - 0.936 g/mL
Acid value : ≤ 4.5
Iodine value : 160 - 200
Saponification value : 188 - 195

Risk and Safety Information

WGK : nwg
RTECS : OI9690000
Flash Point(°F) : 235.4 °F
Flash Point(°C) : 113 °C
Storage Temperature(°C) : Store below 30°C

Transport Information

Marine Pollutant : No
ADR/RID : Not Dangerous Goods
IMDG : Not Dangerous Goods
IATA : Not Dangerous Goods