

RM1300 Glycerol trilaurate, Extra pure**Product Number****Packing**

RM1300 : 1G

Product Information

Product Code : RM1300
Product Name : Glycerol trilaurate, Extra pure
Synonym : Trilaurin, Extra pure; Tridodecanoin
Molecular Formula : $C_{39}H_{74}O_6$
Molecular Weight : 639.00
CAS No. : 538-24-9
EC No. : 208-687-0
HS Code : 2915 90 00
Shelf Life : 4 years

Technical Specification

Appearance : White powder
Solubility : 33.3 mg soluble in 1 mL of chloroform
FTIR (KBr disc) : Matches with the standard pattern
Melting range : ~46.5°C
Assay (HPLC) : min. 99.00 %

Risk and Safety Information

S-Phrase(s) : 22-24/25
WGK : 3
Storage Temperature(°C) : Store at 2 - 8°C

Transport Information

Marine Pollutant : No
ADR/RID : Not Dangerous Goods
IMDG : Not Dangerous Goods
IATA : Not Dangerous Goods