

RM290 Barium chloride dihydrate, Extra pure**Product Number Packing**

RM290	: 500G
RM290	: 5KG


Product Information

Product Code	: RM290
Product Name	: Barium chloride dihydrate, Extra pure
Molecular Formula	: $B a C l_2 \cdot 2 H_2O$
Molecular Weight	: 244.26
CAS No.	: 10326-27-9
EC No.	: 233-788-1
HS Code	: 2827 39 85
EC Index No.	: 056-004-00-8
Shelf Life	: 4 years


Technical Specification

Appearance	: White to off-white crystals or powder
Solubility	: 33.3 mg soluble in 1 mL of water
Calcium (Ca)	: $\leq 0.01\%$
Heavy metals (as Pb)	: $\leq 0.0005\%$
Iron (Fe)	: $\leq 0.0002\%$
Assay (AT/EDTA Titration)	: min. 99.00 %

GHS Safety Information

Hazard Statement(s)	: H301-H332
Precautionary Statement(s)	: P301+P310
Signal Word	: Danger
Hazard Pictogram(s)	

Risk and Safety Information

R-Phrase(s)	: 20-25
S-Phrase(s)	: 45
WGK	: 3
RTECS	: CQ 8751000
Storage Temperature(°C)	: Store below 30°C
Hazard Symbol(s)	

Transport Information

UN No.	: 1564
--------	--------



Product Information

Class	:	6.1
Packaging Group	:	3
Marine Pollutant	:	No
ADR/RID	:	1564 6.1/PG 3
IMDG	:	1564 6.1/PG 3
IATA	:	1564 6.1/PG 3

A-516, Swastik Disha Business Park, Via Vadhani Indl. Est., LBS Marg, Mumbai-400 086, India. Phone:+91-022-6147 1919
Fax:+91-022-6147 1920 Gram : STERILITY . email : info@himedialabs.com. Visit us at our website: www.himedialabs.com

The information contained herein is believed to be accurate and complete. However no warranty or guarantee whatsoever is made or is to be implied with respect to such Information or with respect to any product, method or apparatus referred to herein