

---

## NADP(TPN) disodium salt

**RM392**

---

### Product Identifier

CAS No.	:	24292-60-2
EC No.	:	246-129-8
Molecular Formula	:	$C_{21}H_{26}N_7Na_2O_{17}P_3$
Molecular Weight	:	787.37
Synonym	:	Nicotinamide adenine dinucleotide phosphate disodium salt
HS Code	:	2934 99 00
Storage	:	2 - 8°C
Shelf life	:	3 years

---

### Technical Specification

Appearance	:	White to off-white powder or solid
Solubility	:	33.3 mg soluble in 1 mL of water
Assay (HPLC)	:	min. 90.00%

---

### Safety Information

UN No.	:	Not dangerous goods
Class	:	-
Packing Group	:	-
RTECS	:	Not available
WGK	:	3