
Coenzyme Q10, Hi-LR®

RM6818

Product Identifier

CAS No.	:	303-98-0
EC No.	:	206-147-9
Molecular Formula	:	C ₅₉ H ₉₀ O ₄
Molecular Weight	:	863.34
HS Code	:	2914 62 00
Storage	:	2 - 8°C
Shelf life	:	3 years

Technical Specification

Appearance	:	Faint yellow to dark orange crystals or powder
Solubility	:	33.3 mg soluble in 1 mL of chloroform
Melting range	:	47 - 52°C
Assay (HPLC)	:	>= 98.00%

Safety Information

UN No.	:	Not dangerous goods
Class	:	-
Packing Group	:	-
RTECS	:	DK3900000
WGK	:	1