

Glucoamylase, From *Aspergillus niger*

RM7062

Product Identifier

CAS No.	:	9032-08-0
EC No.	:	232-877-2
EC Index No	:	647-016-00-X
HS Code	:	3507 90 99
Storage	:	2 - 8°C
Shelf life	:	3 years

Technical Specification

Appearance	:	White powder
Solubility	:	Soluble in water.
Activity	:	≥ 30 U/mg solid (salt free)(Starch; pH 4.8; 60°C)

Safety Information

Hazard Pictogram(s)	:	
Signal Word	:	Danger
Hazard Statement(s)	:	H334
Precautionary Statement(s)	:	P261- P284- P304 + P340- P342 + P311- P501
UN No.	:	Not dangerous goods
Class	:	-
Packing Group	:	-
RTECS	:	Not available
WGK	:	1