

Ponceau 2R, Hi-Cert™

RM8826

Product Identifier

CAS No.	:	3761-53-3
EC No.	:	223-178-3
C.I. No.	:	16150
Molecular Formula	:	$C_{18}H_{14}N_2Na_2O_7S_2$
Molecular Weight	:	480.42
HS Code	:	3204 90 00
Storage	:	Below 30°C
Shelf life	:	4 years

Technical Specification

Appearance	:	Red to bordeaux crystals or powder or solid
Solubility	:	20 ppm in water yields clear solution
Absorption maxima	:	501 - 511 nm (20 ppm in water)
Loss on drying (at 105°C, 2 hr)	:	<= 13.00%

Safety Information

Hazard Pictogram(s)



Signal Word	:	Warning
Hazard Statement(s)	:	H351
Precautionary Statement(s)	:	P281
UN No.	:	Not dangerous goods
Class	:	-
Packing Group	:	-
RTECS	:	QJ6825000
WGK	:	3