

Sodium acetate anhydrous, Cell Culture Tested

TC023

Product Identifier

| | | |
|-------------------|---|---|
| CAS No. | : | 127-09-3 |
| EC No. | : | 204-823-8 |
| Molecular Formula | : | C ₂ H ₃ O ₂ Na |
| Molecular Weight | : | 82.03 |
| Synonym | : | Acetic acid sodium salt |
| HS Code | : | 2915 29 90 |
| Storage | : | Below 30°C |
| Shelf life | : | 4 years |

Technical Specification

| | | |
|---------------------------------|---|--|
| Appearance | : | White to off-white hygroscopic crystals or powder or solid |
| Solubility | : | 100 mg soluble in 1 mL of water |
| pH (5% in water at 25°C) | : | 7.00 - 9.20 |
| Loss on drying (at 120°C, 2 hr) | : | <= 1.0% |
| Assay (NT, on dry basis) | : | >= 99.00% |
| Cell Culture Test | : | Passes |

Safety Information

| | | |
|---------------|---|---------------------|
| UN No. | : | Not dangerous goods |
| Class | : | - |
| Packing Group | : | - |
| RTECS | : | AJ4300010 |
| WGK | : | 1 |