



Catalase

Cell Culture Tested

Product Code: TC037

Product Description:

Molecular Weight: 250 kDa

CAS No.: 9001-05-2

Synonym: H₂O₂:H₂O₂ oxidoreductase

Source: Bovine Liver

Quality Control:

Appearance

Brown to dark green solid or powder.

Protein (N x 6.25)

60.00 - 90.00%

Enzymatic Activity units/mg protein

2000.00 - 5000.00 U/mg

Cell Culture Test

Passes

Storage and Shelf Life:

Store at -30°C to -10°C away from bright light.

When frozen, catalase solutions loses the activity, so they should not be frozen or lyophilized.

Direct exposure or solutions to air should be avoided.

Shelf life is 24 months.

Use before expiry date given on the product label.

Disclaimer :

User must ensure suitability of the product(s) in their application prior to use. Products conform solely to the information contained in this and other related HiMedia™ Publications. The information contained in this publication is based on our research and development work and is to the best of our knowledge true and accurate. HiMedia™ Laboratories Pvt Ltd reserves the right to make changes to specifications and information related to the products at any time. Products are not intended for human or animal diagnostic or therapeutic use but for laboratory, research or further manufacturing use only, unless otherwise specified. Statements contained herein should not be considered as a warranty of any kind, expressed or implied, and no liability is accepted for infringement of any patents.

Revision: 03/2025