



## Product Information

### TC040 Chondroitin sulfate A sodium salt, Cell Culture Tested

#### Product Number

#### Packing

TC040	: 25G
TC040	: 5G

#### Product Information

Product Code	: TC040
Product Name	: Chondroitin sulfate A sodium salt, Cell Culture Tested
CAS No.	: 39455-18-0
HS Code	: 3913 90 00
Shelf Life	: 3 years

#### Technical Specification

Appearance	: White to off-white powder
Solubility	: Clear to hazy colourless solution at 10gm in 100ml of water
pH (1% in water at 25°C)	: 5.00 - 7.00
Cell Culture Test	: Passes
Assay (HPLC)	: min. 90.00%

#### Risk and Safety Information

WGK	: 3
Storage Temperature(°C)	: Store at 2 - 8°C

#### Transport Information

Marine Pollutant	: No
ADR/RID	: Not Dangerous Goods
IMDG	: Not Dangerous Goods
IATA	: Not Dangerous Goods

A-516, Swastik Disha Business Park, Via Vadhani Indl. Est., LBS Marg, Mumbai-400 086, India. Phone:+91-022-6147 1919  
Fax:+91-022-6147 1920 Gram : STERILITY . email : [info@himedialabs.com](mailto:info@himedialabs.com). Visit us at our website: [www.himedialabs.com](http://www.himedialabs.com)

The information contained herein is believed to be accurate and complete. However no warranty or guarantee whatsoever is made or is to be implied with respect to such Information or with respect to any product, method or apparatus referred to herein