

**TC041                      Gelatin, Cell Culture Tested****Product Number****Packing**

TC041	: 100G
TC041	: 500G

**Product Information**

Product Code	: TC041
Product Name	: Gelatin, Cell Culture Tested
CAS No.	: 9000-70-8
EC No.	: 232-554-6
HS Code	: 3503 00 10
Shelf Life	: 3 years

**Technical Specification**

Appearance	: White to yellowish brown powder
Solubility	: Hazy faint yellow solution at 5 gms in 100 ml of water, gives clear faint yellow coloured solution on warming
Loss on drying (at 105°C, 2 hr)	: <= 15.00%
Total Nitrogen	: <= 15.50%

**Risk and Safety Information**

WGK	: 3
RTECS	: LX8580000
Storage Temperature(°C)	: Store below 30°C

**Transport Information**

Marine Pollutant	: No
ADR/RID	: Not Dangerous Goods
IMDG	: Not Dangerous Goods
IATA	: Not Dangerous Goods