

Adenosine 5'-triphosphate disodium salt Cell Culture Tested

TC085

Product Identifier

CAS No.	:	987-65-5
EC No.	:	213-579-1
Molecular Formula	:	$C_{10}H_{14}N_5Na_2O_{13}P_3$
Molecular Weight	:	551.14
HS Code	:	2934 99 00
Storage	:	2 - 8°C
Shelf life	:	3 years

Technical Specification

Appearance	:	White to off-white crystals or solid or powder
Solubility	:	33.3 mg soluble in 1 mL of water
FTIR	:	Matches with the standard pattern
Assay (HPLC)	:	min. 95.00%
Cell Culture Test	:	Passes

Safety Information

UN No.	:	Not dangerous goods
Class	:	-
Packing Group	:	-
RTECS	:	Not available
WGK	:	3