

2'-Deoxyadenosine monohydrate, Cell Culture Tested

TC111

Product Identifier

CAS No.	:	16373-93-6
EC No.	:	213-488-7
Molecular Formula	:	C ₁₀ H ₁₃ N ₅ O ₃ .H ₂ O
Molecular Weight	:	269.26
HS Code	:	2934 99 90
Storage	:	2 - 8°C
Shelf life	:	4 years

Technical Specification

Appearance	:	White to beige crystals or powder or solid
Solubility	:	33.3 mg soluble in 1 mL of hot water
Water (K.F.)	:	5.0 - 8.0%
Assay (HPLC/NT, anhydrous basis)	:	99.00 - 102.00%

Safety Information

Hazard Pictogram(s)



Signal Word	:	Danger
Hazard Statement(s)	:	H301
Precautionary Statement(s)	:	P264- P270- P301 + 310- P330- P405- P501
UN No.	:	2811
Class	:	6.1
Packing Group	:	III
RTECS	:	AU7358600
WGK	:	1