
Heparin sodium salt, Cell Culture Tested

TC138

Product Identifier

CAS No.	:	9041-08-1
HS Code	:	3001 90 91
Storage	:	2 - 8°C
Shelf life	:	3 years

Technical Specification

Appearance	:	White to grayish brown amorphous powder
Solubility	:	33.3 mg soluble in 1 mL of water
pH (1% in water at 25°C)	:	5.00 - 7.50
Potency	:	min. 140 USP Units/mg
Cell Culture Test	:	Passes

Safety Information

UN No.	:	Not dangerous goods
Class	:	-
Packing Group	:	-
RTECS	:	MI0850000
WGK	:	2