

## Nicotinic acid

**TC157**

### Cell Culture Tested

### Product Identifier

CAS No.	:	59-67-6
EC No.	:	200-441-0
Molecular Formula	:	C <sub>6</sub> H <sub>5</sub> NO <sub>2</sub>
Molecular Weight	:	123.11
Synonym	:	Niacin; Vitamin B <sub>3</sub>
HS Code	:	2936 29 20
Storage	:	Below 30°C
Shelf life	:	3 years

### Technical Specification

Appearance	:	White to off-white crystals or powder
Solubility	:	10 mg soluble in 1 mL of water
FTIR	:	Matches with the standard pattern
Melting range	:	234 - 240°C
Assay (NaOH Titration )	:	min. 99.00%
Cell Culture Test	:	Passes

### Safety Information

Hazard Pictogram(s)



Signal Word	:	Warning
Hazard Statement(s)	:	H319
Precautionary Statement(s)	:	P305+P351+P338
UN No.	:	Not dangerous goods
Class	:	-
Packing Group	:	-
RTECS	:	QT0525000
WGK	:	1