

---

## Bovine serum Albumin Fraction-V

**TC194**

### Cell Culture Tested

### Product Identifier

---

CAS No.	:	9048-46-8
EC No.	:	232-936-2
HS Code	:	3502 90 00
Storage	:	2 - 8°C
Shelf life	:	3 years

### Technical Specification

---

Appearance	:	White to light brown solid or flakes or powder
Solubility	:	Clear to slightly hazy faint yellow to yellow green solution at 4gm in 100ml of water
pH (1% in water at 25°C)	:	6.50 - 7.50
Loss on drying (at 105°C, 2 hr)	:	<= 5.00%
Assay (Protein, on dry basis)	:	98.00 - 102.00%
Cell Culture Test	:	Complies

### Safety Information

---

UN No.	:	Not dangerous goods
Class	:	-
Packing Group	:	-
RTECS	:	Not available
WGK	:	3