

## Vancomycin hydrochloride

**TC202**

### Cell Culture Tested


#### Product Identifier

CAS No.	:	1404-93-9
Molecular Formula	:	$C_{66}H_{75}Cl_2N_9O_{24}.HCl$
Molecular Weight	:	1485.71
HS Code	:	2941 90 90
Storage	:	2 - 8°C
Shelf life	:	3 years

#### Technical Specification

Appearance	:	White to brown crystals or powder
Solubility	:	100 mg soluble in 1 mL of water
pH (5% in water at 25°C)	:	2.50 - 4.50
Potency	:	$\geq 900 \mu\text{g/mg}$ Vancomycin
Assay (HPLC)	:	min. 80.00%
Cell Culture Test	:	Passes

#### Safety Information

Hazard Pictogram(s)	:	
Signal Word	:	Danger
Hazard Statement(s)	:	H317- H334
Precautionary Statement(s)	:	P261- P272- P280- P284- P302 + P352- P304+P340+P312
UN No.	:	Not dangerous goods
Class	:	-
Packing Group	:	-
RTECS	:	YW4380000
WGK	:	2