



# Sodium bicarbonate

Meets USP-NF, EP, BP, JP and IP testing specifications

### Product Code: TC230M

#### **Product Description :**

Molecular weight: 84.01 Molecular formula: NaHCO<sub>3</sub> CAS No.: 144-55-8

## **Quality Control:**

**Appearance (USP)** White crystalline powder

Appearance (BP, EP) White or almost white crystalline powder

**Appearance (JP)** White crystals or crystalline powder

**Appearance (IP)** A white crystalline powder or small, opaque, monoclinic crystals

**Solubility (USP)** Soluble in water, insoluble in alcohol

**Solubility (BP, EP)** Soluble in water, Practically insoluble in ethanol (96%)

**Solubility (IP)** Freely soluble in water, practically insoluble in ethanol (95%)

**Solubility (JP)** It is soluble in water, and practically insoluble in ethanol (95%) and in diethyl ether.

**pH (5% in water at 25°C) (JP)** 7.90 - 8.40

Identification (JP)

Reactions for sodium salt and for bicarbonate

**Identification A (USP)** Reactions for sodium

**Identification A (BP, EP, IP)** A gas is evolved and solution becomes red

**Identification B (USP, IP)** Reactions for bicarbonate

**Identification B (BP, EP)** Reactions for carbonate and bicarbonate

**Identification C (BP, EP, IP)** Reactions for sodium salts **Appearance of solution (BP, EP, IP)** 5% solution in water is clear and colourless.

**Clarity and colour of solution (JP)** 5% solution in water is clear and colourless.

Insoluble substances (USP) The resulting solution is complete and clear Limit of organics (USP) <= 0.01% Arsenic (USP, EP, BP, JP, IP) <= 0.0002% Carbonate (USP) <= 0.23% Carbonate (EP, BP, IP) The pH of freshly prepared 5% solution in water is not greater than 8.6 Carbonate (JP)

No red colour develops immediately **Normal carbonate (USP)** The solution does not assume more than a faint pink colour immediately

**Chloride (USP)** <= 0.015%

**Chloride (EP, BP)** <= 0.015% **Chloride (JP)** <= 0.040%

**Chloride (IP)** <= 0.02%

**Heavy metals (JP, IP)** <= 0.0005%

Aluminum (USP) <= 0.0002% Iron (USP) <= 0.0005% Iron (EP, BP, IP) <= 0.002% Sulphate (EP, BP) <= 0.015% Sulphate (IP) <= 0.015% Sulphur compounds (USP) <= 0.015% Ammonia (USP) <= 0.002% Ammonium (EP, BP) <= 0.002% Ammonium (JP) The gas evolved does not change moistened red litmus paper to blue Calcium (IP) 2% W/V solution, when boiled for 5 min. is clear Calcium (USP, EP, BP) <= 0.01% Copper (USP) <= 0.0001% Magnesium (USP) <= 0.004% Loss on drying (USP) <= 0.25% (Over silica gel for 4 hr.) Assay (HCl T, dried basis) (USP) 99.00 - 100.50% Assay (HCl titration) (IP, EP, BP) 99.00 - 101.00% Assay (H<sub>2</sub>SO<sub>4</sub> Titration) (JP) min. 99.00%

#### **Storage and Shelf Life:**

Store below 30°C away from bright light. Shelf life is 48 months. Use before expiry date given on the product label.

Disclaimer :

Revision: 04/2023

User must ensure suitability of the product(s) in their application prior to use. Products conform solely to the information contained in this and other related HiMedia<sup>™</sup> Publications. The information contained in this publication is based on our research and development work and is to the best of our knowledge true and accurate. HiMedia<sup>™</sup> Laboratories Pvt Ltd reserves the right to make changes to specifications and information related to the products at any time. Products are not intended for human or animal diagnostic or therapeutic use but for laboratory, research or further manufacturing use only, unless otherwise specified. Statements contained herein should not be considered as a warranty of any kind, expressed or implied, and no liability is accepted for infringement of any patents.