

## DL-Threonine Cell Culture Tested

**TC378**

### Product Identifier

CAS No.	:	80-68-2
EC No.	:	201-300-6
Molecular Formula	:	C <sub>4</sub> H <sub>9</sub> NO <sub>3</sub>
Molecular Weight	:	119.12
Synonym	:	(±)-2-Amino-3-hydroxybutyric acid
HS Code	:	2922 50 90
Storage	:	Below 30°C
Shelf life	:	TC378

### Technical Specification

Appearance	:	White to off-white crystals or powder.
Solubility	:	33.3 mg soluble in 1 mL of water
Melting Point	:	~ 244°C (dec)
Assay (NT)	:	min. 99.00%
Cell Culture Test	:	Passes

### Safety Information

UN No.	:	Not dangerous goods
Class	:	-
Packing Group	:	-
RTECS	:	Not available
WGK	:	1