



Protease

Cell Culture Tested

Product Code: TC401

Product Description:

CAS No. : 9036-06-0

Source : Streptomyces griseus

Storage and Shelf Life:

Store at 2 - 8°C.

Shelf life is 24 months.

Use before the expiry given on the product label.

Quality Control:

Appearance

White to light brown powder.

Solubility

Clear colorless solution at 33.3 mg in 1ml of 0.01M

Sodium acetate with 0.005 M Calcium acetate

at pH 7.5.

Activity

NMT 3.5 Units/mg

Cell Culture Test

Passes

Disclaimer :

Revision: 03/2025

User must ensure suitability of the product(s) in their application prior to use. Products conform solely to the information contained in this and other related HiMedia™ Publications. The information contained in this publication is based on our research and development work and is to the best of our knowledge true and accurate. HiMedia™ Laboratories Pvt Ltd reserves the right to make changes to specifications and information related to the products at any time. Products are not intended for human or animal diagnostic or therapeutic use but for laboratory, research or further manufacturing use only, unless otherwise specified. Statements contained herein should not be considered as a warranty of any kind, expressed or implied, and no liability is accepted for infringement of any patents.