

Forskolin

TC472

Cell Culture Tested

Product Identifier

CAS No.	:	66575-29-9
EC No.	:	266-410-9
Molecular Formula	:	C ₂₂ H ₃₄ O ₇
Molecular Weight	:	410.5
HS Code	:	2932 99 90
Storage	:	at -20°C
Shelf life	:	4 years

Technical Specification

Appearance	:	White to off-white crystals or powder
Solubility	:	33.3 mg soluble in 1 mL of ethanol
Melting range	:	228 - 230°C
Assay (HPLC)	:	min. 99.50%

Safety Information

Hazard Pictogram(s)



Signal Word	:	Warning
Hazard Statement(s)	:	H312
Precautionary Statement(s)	:	P280
UN No.	:	Not dangerous goods
Class	:	-
Packing Group	:	-
RTECS	:	QL6150000