

β-Glucosidase

TC681

From Sweet almond

Cell Culture Tested


Product Identifier

CAS No.	:	9001-22-3
EC No.	:	232-589-7
Molecular Weight	:	~1,10,000
HS Code	:	3507 90 99
Storage	:	2 - 8°C
Shelf life	:	3 years

Technical Specification

Appearance	:	White to tan lyophilized powder
Solubility	:	1.9 - 2.1 mg/mL of water (clear to slight hazy colourless to faintyellow solution)
Activity	:	>=10 U/mg solid
alpha-Amylase	:	<= 0.0005%
Cell Culture Test	:	Passes
Unit definition	:	One unit will liberate 1.0 microgram of glucose per minute at pH 5.0, at 35°C

Safety Information

Hazard Pictogram(s)	:	
Signal Word	:	Danger
Hazard Statement(s)	:	H334
Precautionary Statement(s)	:	P261- P342 + P311
UN No.	:	Not dangerous goods
Class	:	-
Packing Group	:	-
RTECS	:	No available
WGK	:	1