
Glucose-6-Phosphate dehydrogenase

TC682

Cell Culture Tested

Product Identifier

CAS No.	:	9001-40-5
EC No.	:	232-602-6
HS Code	:	3507 90 99
Storage	:	2 - 8°C
Shelf life	:	3 years

Technical Specification

Appearance	:	White amorphous, lyophilized powder.
Solubility	:	Dissolve in 5 mM glycine buffer, pH (with or without 0.1% BSA) or in deionized water at 1 mg/ml (buffer is preferable).
Activity	:	>=400 U/mg solid (D- Glucose-6-phosphate; pH 7.8; 30°C)
Cell Culture Test	:	Passes

Safety Information

UN No.	:	Not Dangerous goods
Class	:	-
Packing Group	:	-
RTECS	:	Not available