
Hexokinase, Cell Culture Tested

TC683

Product Identifier

CAS No.	:	9001-51-8
EC No.	:	232-611-5
HS Code	:	3507 90 99
Storage	:	2 - 8°C
Shelf life	:	3 years

Technical Specification

Appearance	:	White lyophilized powder
Solubility	:	33.3 mg soluble in 1 mL of citrate buffer, pH 7.0.
Activity	:	≥ 140 U/mg protein (Glucose, pH 7.6 at 25°C)
Cell Culture Test	:	Passes

Safety Information

UN No.	:	Not dangerous goods
Class	:	-
Packing Group	:	-
RTECS	:	Not available
WGK	:	3