

Sodium azide

TC704

Cell Culture Tested

Product Identifier

CAS No.	:	26628-22-8
EC No.	:	247-852-1
Molecular Formula	:	NaN ₃
Molecular Weight	:	65.01
HS Code	:	2850 00 30
Storage	:	Below 30°C
Shelf life	:	4 years

Technical Specification

Appearance	:	Colourless or white to off-white solid or crystals or chunks or powder or granules
Solubility	:	Clear colorless solution at 100 mg in 1 mL of water.
Assay (FAS Titration/NT)	:	99.00 - 105.00%
Cell Culture Test	:	Passes

Safety Information

Hazard Pictogram(s)	:	
Signal Word	:	Danger
Hazard Statement(s)	:	H300+H310- H373- H410
Precautionary Statement(s)	:	P273- P280- P301+P310+P330- P302+P352+P310- P391- P501
UN No.	:	1687
Class	:	6.1
Packing Group	:	II
RTECS	:	VY8050000