



Sudan Black B

Cell Culture Tested

Product Code: TC709

Product Description:

Molecular formula: $C_{29}H_{24}N_6$

MW : 456.54

CAS : 4197-25-5

Quality Control:

Appearance

Dark brown to very dark brown to black crystals or powder or solid.

Solubility

20 ppm in ethanol yields clear solution.

Absorption maxima

590.00 - 610.00 (20ppm in ethanol)

Loss on drying (at 105°C, 2 hr)

NMT 5.00%

Cell Culture Test

Passes

Storage and Shelf Life:

Store below 30°C.

Shelf life is 48 months.

Use before the expiry given on the product label.

Disclaimer :

User must ensure suitability of the product(s) in their application prior to use. Products conform solely to the information contained in this and other related HiMedia™ Publications. The information contained in this publication is based on our research and development work and is to the best of our knowledge true and accurate. HiMedia™ Laboratories Pvt Ltd reserves the right to make changes to specifications and information related to the products at any time. Products are not intended for human or animal diagnostic or therapeutic use but for laboratory, research or further manufacturing use only, unless otherwise specified. Statements contained herein should not be considered as a warranty of any kind, expressed or implied, and no liability is accepted for infringement of any patents.

Revision No.: 01/2025