

1-Cyano-4-dimethylaminopyridinium tetrafluoroborate,

TC749

Cell Culture Tested

Product Identifier

CAS No.	:	59016-56-7
Molecular Formula	:	C ₈ H ₁₀ BF ₄ N ₃
Molecular Weight	:	234.99
HS Code	:	2933 39 90
Storage	:	On receipt store at -20°C
Shelf life	:	3 years

Technical Specification

Appearance	:	White to yellow crystals or powder
Solubility	:	33.3 mg soluble in 1 mL of acetonitrile
Melting range	:	194 - 204 °C
Assay (HPLC/Total Nitrogen)	:	min. 97.00%
Cell Culture Test	:	Passes

Safety Information

Hazard Pictogram(s)	:	
Signal Word	:	Warning
Hazard Statement(s)	:	H315- H319- H335
Precautionary Statement(s)	:	P261- P264- P271- P280- P302 + P352- P305+P351+P338
UN No.	:	Not dangerous goods
Class	:	-
Packing Group	:	-
WGK	:	3