

## Krebs-Ringer Bicarbonate Buffer

**With 1.8gms Glucose per litre**

**Without Phenol red, Calcium Chloride and Sodium bicarbonate**

**Product Code: TS1097**

### Product Description:

All media used in tissue culture have a basis of a synthetic mixture of inorganic salts known as a physiological or balanced salt solution (BSS). All the physiological salt solutions have been derived from the salt solution originally described by Sydney Ringer (1885). The first balanced salt solution to be developed specifically for supporting the metabolism of mammalian cells was Tyrode's solution. Since then many modifications have been done to obtain better buffering salt solutions and to prevent calcium precipitation.

The function of a salt solution is:

- To maintain the medium within physiological pH range.
- To maintain intracellular and extra cellular osmotic balance.
- Modified with a carbohydrate, such as glucose serves as an energy source for cell metabolism.

TS1097 is Krebs-Ringer bicarbonate buffer with 1.8gms per litre Glucose. It does not contain phenol red, calcium chloride and sodium bicarbonate.

### Composition:

Ingredients	mg/L
<b>INORGANIC SALTS</b>	
Disodium hydrogen phosphate, anhydrous	100.000
Magnesium chloride anhydrous	46.800
Potassium chloride	340.000
Sodium chloride	7000.000
Sodium dihydrogen phosphate anhydrous	180.000
<b>OTHERS</b>	
D-Glucose	1800.000

### Directions:

1. Suspend 9.5gms in 900ml tissue culture grade water with constant, gentle stirring until the powder is completely dissolved. Do not heat the water.
2. Add 1.260gms of Sodium bicarbonate powder or 16.9ml of 7.5% Sodium bicarbonate solution

(TCL013) for 1 litre of solution and stir until dissolved.

3. Adjust the pH to 0.2-0.3 pH units below the desired pH using 1N HCl or 1N NaOH since the pH tends to rise during filtration.
4. Make up the final volume to 1000ml with tissue culture grade water.
5. Sterilize the solution immediately by filtering through a sterile membrane filter with a porosity of 0.22 micron or less, using positive pressure rather than vacuum to minimize the loss of carbon dioxide.
6. Aseptically dispense the desired amount of sterile solution into sterile containers.
7. Store the liquid solution at ambient temperature and in dark till use.

### Material required but not provided:

Tissue culture grade water (TCL010)  
 Sodium bicarbonate (TC230)  
 Sodium bicarbonate Solution (TCL013)  
 1N Hydrochloric acid (TCL003)  
 1N Sodium hydroxide (TCL002)

### Quality Control:

#### Appearance

Off-white to Creamish white

#### Solubility

Clear solution at 9.5 gms/L

#### pH without Sodium Bicarbonate

6.10 -6.70

#### pH with Sodium Bicarbonate

7.00 -7.60

#### Osmolality without Sodium Bicarbonate (mOsm/Kg H<sub>2</sub>O)

245.00 -285.00

#### Osmolality with Sodium Bicarbonate (mOsm/Kg H<sub>2</sub>O)

255.00 -295.00

**Toxicity test**

Passes

**Endotoxin Content**

NMT 1EU/ml

**Storage and Shelf Life:**

1. All powdered salt mixtures and prepared salt solutions should be stored at ambient temperature. Use before the expiry date. In spite of above recommended storage condition certain powdered salts may show some signs of deterioration / degradation in certain instances. This can be indicated by change in colour, change in appearance and presence of particulate matter and haziness after dissolution.
2. Preparation of concentrated solutions is not recommended as salt complexes having low solubility may precipitate in concentrated solutions.
3. If desired, sterile supplements can be added to the sterile solution observing all sterility precautions. Shelf life of the solution will depend on the nature of supplements added to the solution.

**Disclaimer :**

Revision : 04/2022

User must ensure suitability of the product(s) in their application prior to use. Products conform solely to the information contained in this and other related HiMedia™ publications. The information contained in this publication is based on our research and development work and is to the best of our knowledge true and accurate. HiMedia™ Laboratories Pvt Ltd reserves the right to make changes to specifications and information related to the products at any time. Products are not intended for human or animal or therapeutic use but for laboratory, diagnostic , research or further manufacturing use only, unless otherwise specified. Statements contained herein should not be considered as a warranty of any kind, expressed or implied, and no liability is accepted for infringement of any patents.