
Sulfapyridine

CMS8928

Product Identifier

CAS No.	:	144-83-2
EC No.	:	205-642-7
Molecular Formula	:	$C_{11}H_{11}N_3O_2S$
Molecular Weight	:	249.29
HS Code	:	2918 11 90
Storage	:	Below 30°C
Shelf life	:	3 years

Technical Specification

Appearance	:	White to pale yellow crystals or powder or granules
Solubility	:	Sparingly soluble in acetone
FTIR	:	Matches with the standard pattern
Melting range	:	190 - 193°C
Assay (NaNO ₂ Titration)	:	99.00 - 100.50%

Safety Information

UN No.	:	Not dangerous goods
Class	:	-
Packing Group	:	-
RTECS	:	DA9625000
WGK	:	2