


TO WHOMSOEVER IT MAY CONCERN

Performance Certificate

This is to certify that HiMedia's InstaNX[®] Mag32 System [Cat# LA1096] and MB554MPF32200 HiPurA[®] Multi Sample Pre-filled Plates for Insta NX[®] Mag32 can be used for DNA extraction from different forensic casework samples viz., blood, blood-stained fabric, bone, teeth, nail, tissues, hair shaft, etc. The results demonstrate that DNA extracted on HiMedia's InstaNX[®] Mag32 System [Cat# LA1096] using the MB554MPF32200 HiPurA[®] Multi Sample Pre-filled Plates for Insta NX[®] Mag32 generates STR profile results that are comparable and at par with other extraction systems currently used in DNA Biology Division at Forensic Science Laboratory, Kolkata.

Warm Regards



Assistant Director,
Forensic Science Laboratory, Kolkata,
Govt. of West Bengal

DR. SANDIP GHOSH
Assistant Director
(Biology Division)
Forensic Science Laboratory
Govt. of West Bengal, Kol-37