

Demonstration Report (24th & 25th Aug'2023 at RFSL-Surat)

This is to certify that Demonstration of HiMedia's Instrument and Kits:

- LA1096 - InstaNX®Mag32 (Magnetic Bead Based Extraction System)
- MB554PC16200 - HiPurA®Multi-sample Pre-filled Cartridges for DNA Extraction
- MBPCR141- Hi-PCR® Cow-Buffer Probe PCR Kit
- LA1073 - InstaQ96®Plus (Real-Time PCR)

Was conducted on 24th & 25th Aug'2023 at **Regional Forensics Laboratory (RFSL)-Surat.**

Below Samples types were processed for Extraction using InstaNX®Mag32:

Sr. No.	Sample Type	Source
1	Bone	Human
2	Teeth	Human
3	Teeth	Human
4	Teeth	Human
5	Foetus- Tissue (Formaline)	Human
6	Foetus- Tissue (Fresh)	Human
7	Pant (Fabric)- Blood Stained	Human
8	Pyjama (Fabric)- Blood & Semen Stained	Human
9	Cloth- Meat	Cattle
10	Cloth- Meat	Cattle
11	Cloth- Meat	Cattle
12	Meat (Tissue)	Cattle

Below Samples were processed for Detection using InstaQ96®Plus & MBPCR141:

Sr. No.	Sample Type	Source
1	Cloth- Meat	Cattle
2	Cloth- Meat	Cattle
3	Cloth- Meat	Cattle
4	Meat (Tissue)	Cattle

The result generated for above samples was satisfactory.

(Run File Attached)

Thanking you,
Yours sincerely,

For, HiMedia Laboratories Pvt. Ltd.,

Authorised Signatory



REGISTERED OFFICE -

Plot No. C40, Road No. 21Y, MIDC, Wagle Industrial Area, Thane (West) - 400604, Maharashtra, India.
Tel: +91-22-6147 1919 / 6116 9797 / 6903 4800 | Fax: +91-22-6147 1920
Email: info@himedialabs.com | Web: www.himedialabs.com



... expect only quality from us™
CIN: U85195MH1982PTC028194